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规格: 10mL/100mL/500mL

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Cat: G1061

Size: 10mL/100mL/500mL

Storage: 2-8°C, avoid light, valid for 1 year.

Introduction

Crystal Violet Stain Solution, 2.5% is mainly composed of ammonium oxalate, crystal violet, ethanol etc. The solution is a component of Gram stain. The primary stain is crystal violet, a basic dye taken up by all bacteria due to its ability to rapidly permeate the cell wall. It stains the protoplast of bacteria purple. The potassium-iodine mixture is the mordant which complexes with the primary stain in the cell. In gram-positive cells, the crystal violet-iodine complex is trapped in the cell due to a decrease in cell wall permeability caused by alcohol dehydration. Gram-positive cells are stained purple. In gram-negative cells, the complex is removed by the decolorizer due to an increase in permeability caused by solubility of the lipids in alcohol. The counterstain used must be a contrasting color to the primary stain and, in this case, it is safranin which stains gram-negative cells red.

Protocol(for reference only)

1. Prepare a thin smear on clear and dry glass slide.
2. Air dry and fix it over a gentle flame, while moving the slide in a circular fashion to avoid localized overheating.
3. Flood with Gram's Crystal Violet Solution, 2.5% for 1 min. Wash with tap water.
4. Flood the smear with Gram's Iodine Solution for 1 min.
5. Pour off the Gram's Iodine Solution and gently wash with tap water. Shake off the excess water from the surface.
6. Decolorize with Gram's Decolorizing Solution for 20 to 60 s until the blue dye no longer flows from the smear. Further delay will cause excess decolorization in the gram-positive cells, and the purpose of staining will be defeated.
7. Wash with tap water.
8. Counterstain with Gram's Safranin Solution for 1 min. Wash with tap water.
9. Allow the slide to air dry or blot dry between sheets of clean bibulous paper and view under oil immersion lens.

Result

Gram-Positive Organisms	Bluish Purple
Gram-Negative Organisms	Pinkish Red

Note

1. Before smearing, circle marks should be made on the back in order to judge the position of subsequent tests.
2. The culture time of the bacteria to be tested will affect the staining, and the positive bacteria are often negative when the culture time is too long or they have died or the bacteria are dissolved.
3. For your safety and health, please wear laboratory clothes and disposable gloves.