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Oil Red O Stain Kit, For Cultured Cells

Cat: G1262

Size: 4×20mL/4×50mL

Storage: 2-8°C, avoid light, valid for 1 year.

Kit Components

Reagent	4×20mL	4×50mL	Storage	
Reagent(A): ORO Fixative	20mL	50 mL	RT, avoid light	
Reagent(B): ORO Stain	B1: ORO Stain A	12 mL	30 mL	2-8°C, avoid light
	B2: ORO Stain B	8 mL	20 mL	RT
Mix B1 with B2 as the ratio of 3:2 and place for 10 min to form ORO Stain. It is not suitable to prepare. If there is precipitate filter it before use.				
Reagent(C): Mayer Hematoxylin Solution	20mL	50 mL	2-8°C, avoid light	
Reagent(D): ORO Buffer	20mL	50 mL	RT	

Introduction

Neutral fat stains often use Sudan II, Sudan III, Sudan IV, Sudan black B, Oil red O and so on. Sudan dye is often used in traditional methods. Recently, azo dye oil red O is more suitable for dyeing fat. Oil Red O is a strong lipid solvent and dye. It is easy to bind to triglycerides in droplet shape, but has a slightly poor binding to phospholipids. The dyeing principle is generally regarded as a physical miscibility or adsorption, and the fat is dyed by miscibility. The solubility of dyes in frozen sections is higher than that in the original solvent, so when dyeing, the dyes are transferred from organic solvents to lipids and the fats are dyed.

Oil Red O Stain Kit, For Cultured Cells is shortened as ORO Stain and can show least fatty droplet. The sample should not be fixed with ethanol (10% formalin if required). The positive staining result of fat is orange to red, but the specific color depends on the concentration of lipid.

Self Provided Materials

60% isopropanol, Distilled water

Protocol(for reference only)

For Cultured Cells

1. Remove the medium and gently wash the cells twice with PBS. Add ORO Fixative to the cells and fix for 20-30min.
2. Discard ORO Fixative and wash the cells twice with distilled water.
3. Soak in 60% isopropanol for 5min.
4. Discard 60% isopropanol and soak in ORO Stain for 10-20min.
5. Discard the ORO Stain and rinse in water 2-5 times until no longer stain.
6. Re-dyeing the nucleus with Mayer Hematoxylin Solution for 1-2min. Discard the dye and wash 2-5 times.
7. Rinse the section in ORO Buffer for 1min and discard.
8. Cover the cells with water and view under the microscope.

For Cell Smear

1. Prepare fresh bone marrow and blood smear, then fix in ORO Fixative for 10-15 min.
2. Take out the smear and put it in the air for 10-15 min.
3. Soak the smear in ORO Stain and dye for 15min.
4. Rinse with 60% isopropanol for 20-30s, rinse with running water, and wash slightly with distilled water. Re-dyeing the nucleus with Mayer Hematoxylin Solution for 1-2 min.
5. Soak in ORO Buffer for 1min.

Result

Neutral Fat	Orange or Red
Phospholipid	Pink Red
Nucleus	Blue

Note

1. ORO Stain is not stable enough and easy to precipitate, so it is not suitable to prepare it in advance.

2. Mayer's Hematoxylin Solution should not stain for too long.
3. The dyeing results can not last for a long time, and should be viewed and photographed as soon as possible.
4. For your safety and health, please wear experimental clothes and disposable gloves.