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规格: 50mL /100mL/500mL

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Cat:G1286

Size:50 mL/100mL/500mL

Storage:2-8℃, avoid light, valid for 6 months.

Introduction

Schiff reagent is often used with periodic acid, the method is called PAS Stain(Periodic Acid-Schiff Stain) and mainly used to detect glycogen or other polysaccharide substances in tissues. Periodic acid is a kind of oxidant, which can oxidize the adjacent diol groups in the polysaccharide molecule to dialdehyde groups. The aldehyde groups can react with Schiff Reagent to form red insoluble complex.

Protocol (for reference only)

1. Fix sections in 10% formalin fixative. Conventionally dehydration and embedding.
2. For paraffin section, dewax to distilled water; For frozen section, restore room temperature.
3. Rinse in tap water for 2-3min, and then wash with distilled water twice.
4. Place in Oxidant for 5-8min at RT, generally not more than 10min.Wash with distilled water twice.
5. Dye with Schiff Reagent in a dark place at RT, and stain for 10-20min. Rinse in tap water for 10min.
6. Stain with Mayer Hematoxylin Solution for 1-2min.
7. Differentiate by Acidic Differentiation Solution for 2-5s. Wash with tap water for 10-15 min to blue.
8. Conventional dehydration by series of ethanol. Transparent by xylene and seal with resinene.

Result

PAS Reaction Positive Substance	Red or Purplish Red
Nucleus	Blue
Cytoplasm	Red in different degrees

Note: the color depth depends on the time that the sample has been in the Oxidant and Schiff Reagent.

Negative Control(optional)

1. Take 1g of amylase and dissolve it in 100ml of PBS (pH5.3), treat it for 30-60 min, and add into Oxidant together with other sections. The result should be negative.
2. (alternative) take the saliva slice (after filtration) and treat it for 30-60 min, then add into Oxidant together with other slices. The result should be negative.

Note

1. For your safety and health, please wear experimental clothes and disposable gloves.
2. Reagents should be stored at low temperature and restore room temperature before use.
3. The incubation temperature of periodate acid should not be higher than 20 °C. When the room temperature is high, shorten the incubation time appropriately. The optimal dyeing temperature of Schiff Reagent is 37℃, adjusting incubation time with low or high temperature.
4. It is recommended to take 1-2 samples for preliminary experiment before formal experiment .

Related Products

G1010 Giemsa Stain Solution(Working Solution)

G1040 Wright Stain Solution

G1120 Hematoxylin-Eosin (HE) Stain Kit

G1080 Mayer 's Hematoxylin Stain Solution(For IHC)

G1260 Oil Red O Saturated Solution, 0.5%, Isopropanol Solvent