## Ehrlich 试剂

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规格: 100mL/500mL

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**Cat:** G1290

**Size:** 100mL/500mL

**Storage:** RT, avoid light, valid for 6 months.

### Introduction

The active ingredient of Ehrlich Reagent is para-Dimethylaminobenzaldehyde. The principle of Ehrlich Reagent is that some bacteria with tryptophanase can decompose tryptophan in peptone to product indole (matrix). Then indole binds to para-Dimethylaminobenzaldehyde to form the red rosindole. Ehrlich Reagent is especially suitable for the identification of Enterobacteriaceae, non-fermenting bacteria, foster bacteria and anaerobic bacteria.

### Self provided materials

Test tube, Peptone Medium, Constant temperature incubator

### Protocol(for reference only)

1. Inoculate the bacteria into the test tube with peptone medium and culture it in a constant temperature incubator at 35-37°C for 24-48 h.
2. Slowly add a small amount of ether or xylene along the tube wall, shake the tube to extract and concentrate indole, then rest until ether or xylene floats on the surface of the medium.
3. Slowly add a few drops of Ehrlich Reagent along the tube wall to observe the color change of the ether or xylene layer.

### Result

Positive	Red
Negative	No change in color

### Note

1. When culturing bacteria, pay attention to avoid pollution.
2. Ehrlich Reagent is only used for scientific research purposes, but not for clinical diagnosis.
3. For your safety and health, please wear experimental clothes and disposable gloves.