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Prussian Blue Iron Stain Kit(With Neutral Red Solution)

Cat: G1420

Size: 2×50mL/2×100mL

Storage: RT, avoid light, valid for 1 year.

Introduction

Hemosiderin is a hemoglobin derived pigment, which is golden yellow or brownish yellow particles. Because it contains iron and golden yellow, it is called hemosiderin. When the red blood cells are engulfed by macrophages, under the action of lysosomal enzymes, hemoglobin is broken down into iron free orange blood and iron-containing hemosiderin.

Perls Prussian blue reaction, also known as hemosiderin staining, can produce blue after being treated with potassium ferrocyanide and dilute acid, which is common in the interstitium of phagocytes, mainly showing ferric iron salts. Perls Prussian blue reaction is a very classical histochemical reaction, which is a sensitive and traditional excellent method to display the ferric iron in tissues. Its dyeing principle is: potassium ferrocyanide solution separates the ferric iron from the protein by dilute hydrochloric acid, and the ferric iron reacts with potassium ferrocyanide to form an insoluble blue compound named Prussian blue. Ferrocyanide of ferric iron is a very stable compound, which can be re-dyed with red dye after reaction, such as nuclear fast red, eosin, neutral red, etc.

Perls Stain is often used to display various hemorrhagic lesions in local tissues, and it is common in phagocytes. Perls reaction can be used to determine the deposition of hemosiderin, and this staining method can distinguish hemosiderin from other pigments. The dyeing solution has good stability, can be preserved for a long time, is not easy to produce precipitation, has a wide range of applications, and can be used for re-dyeing.

Kit Components

Reagent		2×50mL	2×100mL	Storage
Reagent(A): Perls Stain	A1: Perls Stain A	25mL	50mL	RT, avoid light
	A2: Perls Stain B	25mL	50mL	RT
Before use, mix equal parts of A1 and A2 to form Perls Stain. It is not suitable to prepare in advance.				
Reagent(B): Neutral Red Solution		50mL	100mL	RT, avoid light

Self Provided Materials

10% neutral formalin fixative, Series of ethanol, Distilled water, 4% paraformaldehyde

Protocol(for reference only)

(一) For paraffin section staining

1. Fix the tissue in 10% neutral formalin fixative, dehydrate and embed.
2. Cut section in 4um thick and dewax to distilled water. Rinse in distilled water for 1min.
3. Soak the section in Perls Stain(see note 4) and stain for 15-30min. Rinse fully in distilled water for 2-5min.
4. Light staining the nucleus with Neutral Red Solution for 15-30s. Rinse in tap water for 1-5s.
5. Conventionally dehydrate and transparent, seal with resinene.

(二) For frozen section staining

1. Without dewaxing, rinse directly and quickly with distilled water for 2-3min.
2. Follow the other steps as paraffin sections.

(三) For cultured cell staining

1. Fix in 4% paraformaldehyde for 10-20min.
2. Rinse in tap water twice for each time 2min.
3. The steps of Staining, dehydration, transparency and sealing are the same as the steps of paraffin section. The time should be shortened accordingly.

Result

Hemosiderin or Ferric Iron	Blue
Nucleus and other Tissues	Red

Negative Control

Take the same adjacent section and dewax to water. After incubation in 5% oxalic acid for 2-6h, the

procedure is the same as above. The result should be negative.

Note

1. Section dewaxing should be as clean as possible.
2. 10% neutral formalin is often used for tissue fixation. After long-term fixation with common formalin, tissue will be damaged. Avoid the use of acid fixatives, chromate treatment will also hinder the preservation of iron.
3. During the whole operation process, the container should be clean and avoid the use of metal iron products. When washing sections and containers, distilled water is suitable, because ordinary water contains iron.
4. When dyeing with Perls Stain, the time should be adjusted according to the sample situation.
5. All sections should use the same positive control section, so it is very important to select the appropriate control. Autopsy lung tissue is a good control, containing a considerable number of iron positive macrophages (heart failure cells).
6. Series of ethanol should be replaced frequently.
7. For the staining of frozen section and cell, should explore the experimental conditions according to the specific conditions.
8. For your health and safety, please wear the experimental clothes and disposable gloves.