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Prussian Blue Iron Stain Kit (With Eosin Solution)

Cat: G1424

Size: 2×50mL/2×100mL

Storage: RT, avoid light, valid for 1 year.

Introduction

Perls Stain is often used to display various hemorrhagic lesions in local tissues, and it is common in phagocytes. Perls reaction can be used to determine the deposition of hemosiderin, and this staining method can distinguish hemosiderin from other pigments. The dyeing solution has good stability, can be preserved for a long time, is not easy to produce precipitation, has a wide range of applications, and can be used for re-dyeing.

In the Prussian Blue Iron Stain Kit (With Eosin Solution), eosin is used as the counterstain solution. Eosin is a common dye solution, and the dyeing time of this solution is shorter than that of nuclear fast red.

Kit Components

Reagent		2×50mL	2×100mL	Storage
Reagent(A):	A1: Perls Stain A	25mL	50mL	RT, avoid light
Perls Stain	A2: Perls Stain B	25mL	50mL	RT
Before use, mix equal parts of A1 and A2 to form Perls Stain. It is not suitable to prepare in advance.				
Reagent(B):	Eosin Solution	50mL	100mL	RT, avoid light

Protocol(for reference only)

(一) For paraffin section staining

- Dewax to distilled water
Xylene twice each for 5-10min. 100% ethanol twice each for 3-5min.
95%, 85%, 75% ethanol each for 3-5min. Rinse in tap water or distilled water for 1-3min.
- Stain
Stain with Perls Stain(see note 4) for 15-30min. Rinse in distilled water for 5-10min.
Light re-dyeing with Eosin Solution for 15-30s. Rinse quickly in tap water for 2-3s.
- Dehydrate, transparent and seal
75%, 85%, 95% ethanol each for 3-5s.
100% ethanol twice each for 1min.
Xylene twice each for 1-2min. Seal with resinene.

(二) For frozen section staining

- Without dewaxing, rinse directly and quickly with distilled water for 2-3min.
- Follow the other steps as paraffin sections.

(三) For cultured cell staining

- Fix in 4% paraformaldehyde for 10-20min.
- Rinse in distilled water twice for each time 2min.
- The steps of staining, dehydration, transparency and sealing are the same as the steps of paraffin section. The time should be shortened accordingly.

Result

Hemosiderin or Ferric Iron	Blue
Nucleus and other Tissues	Red

Negative Control

Take the same adjacent section and dewax to water. After incubation in 5% oxalic acid for 2-6h, the procedure is the same as above. The result should be negative.

Note

- Section dewaxing should be as clean as possible.
- 10% neutral formalin is often used for tissue fixation. After long-term fixation with common formalin, tissue will be damaged. Avoid the use of acid fixatives, chromate treatment will also hinder the preservation of iron.
- During the whole operation process, the container should be clean and avoid the use of metal iron products. When washing sections and containers, distilled water is suitable, because ordinary water contains iron.

4. When dyeing with Perls Stain, the time should be adjusted according to the sample situation.
5. All sections should use the same positive control section, so it is very important to select the appropriate control. Autopsy lung tissue is a good control, containing a considerable number of iron positive macrophages (heart failure cells).
6. Series of ethanol should be replaced frequently.
7. For the staining of frozen section and cell, should explore the experimental conditions according to the specific conditions.
8. For your health and safety, please wear the experimental clothes and disposable gloves.