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**Cat:** G1511

**Size:** 3×50mL

**Storage:** RT, avoid light, valid for 1 year.

### Kit Components

Reagent	3×50mL	Storage
Reagent(A): Mayer Hematoxylin Staining Solution	50mL	RT, avoid light
Reagent(B): Sudan III Differentiation Solution	50mL	RT
Reagent(C): Sudan III Solution	50mL	RT, avoid light

### Introduction

Neutral fat stains often use Sudan II, Sudan III, Sudan IV, Sudan black B, oil red O and so on. The dyeing principle is generally regarded as a physical miscibility or adsorption, and the fat is dyed by miscibility. The solubility of dyes in frozen sections is higher than that in the original solvent, so when dyeing, the dyes are transferred from organic solvents to lipids and the fats are dyed.

Sudan III Stain Kit is mainly used to show fatty degeneration of tissues and organs and abnormal lipid-like sedation. It often occurs in fatty degeneration of liver, kidney, heart and other parenchymal organs that most of the neutral fat droplets appear in cells. It can identify and diagnose tumors and their properties in adipose tissue. The sample should not fixed with ethanol (10% formalin if required) or paraffin. Frozen sections or carbon wax sections were needed.

### Self Provided Materials

70% ethanol, Distilled Water, Microscope, Glycerol Gelatin.

### Protocol (for reference only)

1. Prepare low-temperature sections of fresh tissue, generally at  $-20^{\circ}\text{C}$  to  $-25^{\circ}\text{C}$ . If the sample is lipoma, it should be adjusted to  $-30^{\circ}\text{C}$ .
2. Slice the frozen section at 6-15  $\mu\text{m}$  (6-8  $\mu\text{m}$  is preferable) and attach to the slide. Wash with distilled water.
3. Re-dyeing with Mayer Hematoxylin Staining Solution for 1 min.
4. After tap water washing, differentiate by Sudan III Differentiation Solution for several seconds, and wash with tap water till the nucleus is blue.
5. Rinse with distilled water and soak slightly with 70% ethanol.
6. Stain with Sudan III Solution for 30min.
7. Differentiate with 70% ethanol for several seconds and wash in tap water.
8. Use filter paper to absorb the slices and surrounding water and let them dry slightly.
9. Glycerol gelatin sealing.

### Result

Neutral fat	Orange Red
Nucleus	Light Blue

### Note

1. Samples should not be fixed with ethanol or embedded with paraffin, and frozen sections should be used. It is recommended to use 10% neutral formalin or 10% formaldehyde-calcium solution to fix.
2. Precipitation of dyes must be prevented during dyeing. Therefore, slices should be covered when they are put into dye solution, so as not to contact with flowing air to avoid precipitation when the solution volatilizes.
3. Frozen sections are easier to be colored, and overdyeing should be avoided when re-dyeing with Mayer hematoxylin.
4. Sudanese dyes fade easily and should be kept in airtight condition.
5. Samples sealed with glycerol gelatin did not last long. If long-term preservation is required, the edge of the boundary between the cover slide and the slide can be sealed with neutral gum.
6. For your safety and health, please wear lab clothes and disposable gloves.