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**Cat:** G1560

**Size:** 3×50mL/3×100mL

**Storage:** RT, avoid light, valid for 6 months.

### Introduction

Alcian blue is a kind of copper titanium cyan conjugate dye, which was originally used for textile fiber dyeing. This kind of cationic dye combines with acid group, that is, alcian blue forms insoluble complex with anionic group such as carboxyl group and sulfate group. Alcian blue is composed of copper containing phthalocyanine ring in the center and four isothioureia groups connected by thioether bond. The isothioureia group is moderately alkaline, which makes alcian blue cationic.

According to the different pH value of dye solution, the genus of mucus can be determined. When pH value is 1, the carboxyl group cannot be ionized and thus cannot dye, but the sulfuric acid group can be shown. When pH value is 2.5, the mucin with carboxyl group (such as proteoglycan / hyaluronic acid and epithelial acid mucin) is well dyed, while the sulfated mucin is poorly dyed. The neutral mucin in gastric mucosa and Brunner gland could not react with alcian blue. It is often used to distinguish mucinous epithelial tumors and to prove whether there is mucinous substance in tumors.

### Kit Components

Reagent	3×50mL	3×100mL	Storage
Reagent (A): Alcian Acidic Solution	50mL	100mL	RT
Reagent (B): Alcian Staining Solution	50mL	100mL	RT, avoid light
Reagent (C): Nuclear Fast Red Solution	50mL	100mL	RT, avoid light

### Protocols(for reference only)

1. For paraffin sections, dewax to distilled water. For cell smear or cell suspension, fix in paraformaldehyde.
2. Soak in Alcian Acidic Solution for 3mins.
3. Stain with Alcian Staining Solution for 30mins.
4. Wash in running tap water.
5. Re-dyeing with Nuclear Fast Red Solution for 5mins.
6. Wash in running tap water for 1 min.
7. Dehydrate in series of alcohol, transparent by xylene and seal with resinene.

### Result

Sulfated Mucin (Sulfated Mucin and Salivary Mucin)	Blue
Carboxylated Mucin (Proteoglycan, Hyaluronic Acid)	Blue
Nucleus	Red

### Note

1. The fixative adopts 10% neutral formalin.
2. To selectively identify sulfated mucin and carboxylated mucin, alcian blue with pH is 1.0 should be used, it is to adjust the pH value to 1.0. The dyeing time can be adjusted according to the dyeing degree.
3. This reagent is only suitable for scientific research, not for clinical diagnosis or other purposes.
4. For your safety and health, please wear experimental clothes and disposable gloves.

### Related Products

G1563 Alcian Blue Stain Kit, pH 1.0

G1562 Alcian Blue Solution, pH 2.5