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规格: 100mL /500mL

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**Cat:** G1615

**Size:** 100mL /500mL

**Storage:** RT, avoid light, valid for 6 months.

### Introduction

Papanicolaou method is widely used in routine cytological staining. Papanicolaou stain is initially used to detect estrogen levels in vaginal epithelium and pathogens such as Candida and trichomonad in the genital tract. The combination of orange G6 and EA36 or EA50 can dye the cytoplasm into bright green, blue and pink. At present, most laboratories use finished dye solution, so we should pay attention to the best conditions of each dye solution after improvement. The final cytoplasmic staining should be transparent, and the nuclear chromatin should be easily distinguished.

At present, the improved pasteurization solution contains many kinds of ions, which has polychromatic dyeing efficiency. After staining, the cytoplasm is bright and transparent, and the structure of nuclear membrane, nucleolus and chromatin is clear. The nuclear staining solution is mainly Harris hematoxylin staining solution, and the cytoplasm staining solution is mainly EA36 staining solution and EA50 staining solution. Papanicolaou Stain is used for exfoliated cells. The nucleus is blue or black, and the cytoplasm of keratinized squamous cells is pink or orange red.

Papanicolaou EA65 Stain Kit use EA65 staining solution for cytoplasmic staining and non-toxic modified hematoxylin staining solution for nuclear staining. It is more suitable for non gynecological cytology smear staining.

### Self Provided Materials

Fixative (like AF Fixative), Series of alcohol, Microscope, Acid alcohol differentiation solution

### Protocol (for reference only)

1. For cell smear, fix with 95% alcohol for 10-15mins.
2. Rinse in 95%, 85%, 75% alcohol separately for 1min. Rinse with tap or distilled water for 1min.
3. Stain with Hematoxylin Solution for 5-10mins. Rinse with tap water for 2mins.
4. Differentiate with 1% acid alcohol differentiation solution for about 4-5s or 0.5% acid alcohol differentiation solution for about 10s. Rinse with tap water for 2mins.
5. Blue with Bluing Solution for 2mins. Rinse with tap water for 2mins.
6. Dehydrate in 75%, 85%, 95% alcohol separately for 2mins.
7. Stain with Orange G Solution for 2mins. Dehydrate in 95% alcohol for 2mins twice.
8. Stain with EA65 Solution for 2mins. Dehydrate in 95% alcohol for 2mins twice.
9. Dehydrate in absolute alcohol for 1min twice.
10. Transparent with xylene and seal with resinene.

### Result

Nucleus	Blue Purple or Black
Cytoplasm of non keratinized cells	Light Blue or Light Green
Cytoplasm of keratinocytes	Pink or Orange

### Note

1. All dye solutions should be filtered before use and changed frequently.
2. For your safety and health, please wear experimental clothes and disposable gloves.