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**Cat:** G1690

**Size:** 2×50mL

**Storage:** RT, avoid light, valid for 6 months.

### Kit Components

Reagent	2×50mL	Storage
Reagent (A): Sudan Black B Stain	50mL	RT, avoid light
Reagent (B): Nuclear Fast Red Stain	50mL	RT, avoid light

### Introduction

Neutral fat stains often use Sudan II, Sudan III, Sudan IV, Sudan black B, oil red O and so on. Sudan is often used in traditional methods. Recently, azo dye oil red O is more suitable for dyeing fat. Oil red O is a strong lipid solvent and dye. It is easy to bind to triglycerides in droplet shape, but has a slightly poor binding to phospholipids.

Sudan Black B Stain Kit is mainly used to show fatty degeneration of tissues and organs and abnormal lipid-like sedation. It often occurs in fatty degeneration of liver, kidney, heart and other parenchymal organs that most of the neutral fat droplets appear in cells. It can identify and diagnose tumors and their properties in adipose tissue. The sample should not fixed with ethanol (10% formalin if required) or paraffin. Frozen sections or carbon wax sections were needed. After Sudan Black B staining, lipoma, liposarcoma, follicular membrane tumor and adrenal cortex adenoma show positive reaction, while renal clear cell carcinoma and ovarian fibroma show negative reaction.

### Self Provided Materials

Slide, 70% ethanol, Distilled water, Microscope, Glycerol gelatin.

### Protocol (for reference only)

1. Prepare low-temperature sections of fresh tissue, generally at  $-20^{\circ}\text{C}$  to  $-25^{\circ}\text{C}$ . If the sample is lipoma, it should be adjusted to  $-30^{\circ}\text{C}$ .
2. Slice the frozen section at 5-10  $\mu\text{m}$  (6-8 micron is preferable) and attach to the slide.
3. Soak slightly with 70% ethanol.
4. Stain with Sudan Black B Stain for 1-2mins. Wash with 70% ethanol to remove excess solution.
5. Rinse with distilled water for 1 min.
6. Re-dyeing slightly with Nuclear Fast Red Stain for 5-10mins.
7. (optional) Under the view of microscope, if the staining is too deep, can differentiate by acid ethanol differentiation solution for several seconds.
8. Rinse with tap water and wash with distilled water.
9. Use filter paper to absorb the slices and surrounding water and let them dry slightly.
10. Seal with glycerol gelatin or arabinose gum.

### Result

Neutral fat	Black
Phospholipid	Dark Gray
Nucleus	Red

### Note

1. Samples should not be fixed with ethanol or embedded with paraffin, and frozen sections should be used. It is recommended to use 10% neutral formalin or 10% formaldehyde-calcium solution to fix.
2. Precipitation of dyes must be prevented during dyeing. Therefore, slices should be covered when they are put into dye solution, so as not to contact with flowing air to avoid precipitation when the solution volatilizes.
3. Sudanese dyes fade easily and should be kept in airtight condition.
4. Samples sealed with glycerol gelatin did not last long. If long-term preservation is required, the edge of the boundary between the cover slide and the slide can be sealed with neutral gum.
5. For your safety and health, please wear lab clothes and disposable gloves.