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规格: 100mL/500mL

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Cat: G1810

Size: 100mL/500mL

Storage: RT, avoid light, valid for 1 year.

Introduction

Evans Blue and trypan blue are both cell reactive dyes, which are commonly used to detect the integrity of cell membrane and cell survival. Living cells can not be dyed blue by Evans Blue because of their efflux function, while dead cells can be dyed light blue. Therefore, dead cells and living cells can be distinguished under microscope by this method, but death and necrosis cannot be distinguished.

Protocols(for reference only)

Blood-brain barrier permeability

1. Take the treated experimental animals (such as mice), after intravenous injection of Evan's Blue Stain Solution, 0.5% , the eyes and skin of mice appear blue. After 0.5-1 h, kill the mice and pick the brain tissues.
2. Place the brain tissue in a 1.5mL centrifuge tube and add 1mL PBS. Quickly make the brain tissue into a homogenate with a tissue homogenizer and centrifuge.
3. Take the supernatant, add equal amount of trichloroacetic acid, and incubate at 4 °C. This step can also be operated as follows: take the supernatant, add acetone according to the ratio of supernatant: acetone is 3:7, and incubate at room temperature for 24h.
4. Centrifugation for 15mins.
5. Take the above solution and measure the absorbance value (OD value) at 620nm with a spectrophotometer. At the same time, measure the OD values of standard Evans Blue with different grades and draw the standard curve. Calculate The Evans Blue content of the sample to be tested according to the standard curve.

Living cells staining

1. Take 100μl of re-suspension cells into the conventional centrifuge tube, add 100μl Evan's Blue Stain Solution, 0.5% into the tube and mix it gently for staining (the staining time can be appropriately prolonged, but not more than 10 minutes).
2. Take a small number of stained cells and count them with a blood cell counting board. Usually, if you want to make a more accurate quantitative analysis, count at least 500 cells for each cell sample and count the number of blue cells and the total number of cells. The cell survival rate is calculated as follows:
Cell survival rate = (total number of cells - number of blue cells) / total number of cells × 100%

Note

1. Evan's Blue Stain Solution, 0.5% is slightly toxic to human body, pay attention to self protection.
2. When staining cells, note that apoptotic bodies occasionally reject staining.
3. In the blood-brain barrier permeability experiment, the injection amount of Evan's Blue Stain Solution, 0.5% should be adjusted according to the weight of different animals and animals.
4. It is best to use a cryocentrifuger for centrifugation.
5. For your safety and health, please wear experimental clothes and disposable gloves.