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Cat: G2110

Size: 2×20mL

Storage: -20℃, avoid light, valid for 6 months.

Kit Components

	Reagent	2×20mL	Storage
Reagent (A): AChE Incubation Solution	A1: AChE Iodide	10mg	4℃, avoid light
	A2: AChE Buffer A	2mL	RT
	A3: AChE Buffer B	14mL	4℃, avoid light
	A4: AChE Buffer C	2mL	4℃, avoid light
	A5: AChE Buffer D	2mL	4℃, avoid light
	A6: Iso-OMPA	0.4mL	-20℃, avoid light
Reagent (B):Hematoxylin Solution		20mL	4℃, avoid light
Reagent (C): AChE-ChE Inhibitor		0.4mL	4℃, avoid light

Introduction

Cholinesterase(ChE) belongs to specific esterase and can be divided into two categories. One is acetylcholinesterase, also known as true cholinesterase, which can hydrolyze acetylcholine and play a physiological regulatory role; the other is cholinesterase, also known as Pseudo cholinesterase(PsChE), which can hydrolyze choline esters instead of acetylcholinesterase. Acetylcholinesterase mainly exists in the cytoplasm of neurons, the junction of nerves and muscles which also called motor endplate; PsChE mainly exists in plasma, pancreas and salivary gland, but its physiological function is not clear.

Acetylcholinesterase Stain Kit(Copper Ferricyanide) belongs to Karnovsky and Roots method. Its dyeing principle is that acetylcholinesterase hydrolyzes acetylcholine iodide and releases acetic acid and thiocholine. The thiohydrogen group (-SH) in thiocholine reduces potassium ferricyanide to potassium ferrocyanide, which combines with copper ions to form insoluble red brown to dark brown copper ferrocyanide precipitated in the active site of the enzyme. Its advantages are simple operation, less esterase diffusion, and its disadvantages are poor permeability of substrate to tissue. The staining solution can be used to observe the changes under central and peripheral nerve fibers diseases, and it is also helpful for the diagnosis of megacolon and intestinal neuron dysplasia. When poisoned by organic pesticides, the enzyme is inhibited and the activity of the enzyme decreased with the negative reaction.

Self Provided Materials

10% Formaldehyde Calcium Fixative, Constant temperature incubator, Optical microscope

Protocol(for reference only)

1. Cut frozen section in 6μm thickness, unfix or fix in precooled 10% Formaldehyde Calcium Fixative for 10mins.
2. Wash with distilled water for 3 times and each time for 1 min.
3. Prepare AChE Incubation Solution: before use, take A2 and add into A1 to make the latter solution completely dissolved to form A12 mixture, and store at 4℃. Take appropriate amount of A12 mixture, A3, A4, A5, A6, and mix them fully as the ratio of 1:7:1:1:0.2 to form AchE Incubation Solution, which shall be used within 6h.

Note: if you want to display AChE and ChE, you do not need to distinguish AChE and ChE, and you do not need to add Iso-OMPA.

4. Take the section into the preheated AChE Incubation Solution and incubate at 37℃ in dark for 1-3 h (generally no more than 6 h), then take out when the section is light brown.
5. Wash with distilled water. View under the optical microscope, if the color of active site is still light, can incubate after washing with distilled water until the reaction is appropriate.
6. Rinse with running water for 5mins.
7. Redyeing with Hematoxylin Solution for 3-5mins.
8. Rinse with running water for 10mins.

9. Conventionally dewax in alcohol and transparent by xylene, then seal with resinene.

Result

Active site of AchE enzyme	Red Brown to dark Brown
Nucleus	Blue

Negative control(optional): Take the prepared AChE Incubation Solution and mix it with AChE-ChE Inhibitor as the ratio of 50:1. Take the same section and put it into the AChE Incubation Solution containing AChE inhibitor, then follow the steps as the same as above. The result is negative reaction.

Note

1. The staining solution is suitable for frozen sections, and the exposure time of sections at room temperature should be reduced.
2. For your safety and health, please wear experimental clothes and disposable gloves.