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## Alkaline Fast Green Stain Solution

**Cat:** G2551

**Size:** 2×50mL

**Storage:** RT, avoid light, valid for 1 year.

### Kit Components

Reagent		2×50mL	Storage
Reagent(A): Acid Differentiation Solution		50mL	RT, avoid light
Reagent(B): Alkaline Fast Green Solution	B1: Alkaline Fast Green Solution A	25mL	RT, avoid light
	B2: Alkaline Fast Green Solution B	25mL	RT, avoid light
Before use, mix B1 with B2 in equal amount to form Alkaline Fast Green Solution. It is ready to use.			

### Introduction

Different amino acids have side chain groups with different chemical properties. Some have basic side chains, some have acid side chains, and the resulting proteins have different numbers of basic and acid groups. These groups will make proteins with different net charges in different pH solutions, and the whole protein molecule has more positive charges, that is, basic protein (isoelectric point tends to be alkaline); The whole protein has more negative charge, that is acid protein (isoelectric point tends to be acid).

The principle of Alkaline Fast Green Stain Solution is that basic protein combines with basic dye containing negative charge to show color. The most abundant basic protein in the cell is histone, which mainly exists in the nucleus. So after dyeing, most of the nucleus are dyed green. The staining solution is only used in scientific research field, and is not suitable for clinical diagnosis or other purposes.

### Self Provided Materials

Slide, 70% ethanol, Water bath, Microscope.

### Protocol (for reference only)

#### For Tissue Staining

1. Fix the tissue in 10% formalin and routinely dehydrate and embed. Slice into 4μm. Conventional dewaxing to water.
2. Immerse the section in Acid Differentiation Solution, soak in 60°C water for 30 min, and wash twice with distilled water for 1 min each time.
3. Filter paper absorbs residual water.
4. Dip the section into Alkaline Fast Green Solution for 5-15 min and wash with distilled water for 30 s.
5. Dry at room temperature and seal with neutral gum.

#### For Adherent Cells or Climbing Section Staining

1. The cells were fixed in neutral formalin or tissue cell solution at room temperature for 10-20 min.
2. Add Acid Differentiation Solution according to half of the added amount of medium, incubate in warm at 60 °C for 30min, and wash twice with distilled water for 1min each time. 3. Absorb excess water.
3. Add the same amount of Alkaline Fast Green Solution to dye for 5-15 min, and wash with distilled water for 30 seconds.
4. Observe with water.

#### For Blood Smear or Suspension Cell Staining

1. Take 1 drop of fresh blood to one end of the slide, push the slide and dry at room temperature.
2. Immerse the smear in 70% ethanol for 5min and dry it at room temperature.
3. Immerse the smear in Acid Differentiation Solution and incubate in 60 °C water bath for 30min.
4. Rinse with running water completely and absorb residual water by filter paper .
5. Immerse the smear in Alkaline Fast Green Solution for 15-30 min.
6. Rinse with running water and dry at room temperature.
7. View under the microscope directly or after sealing with resinene.

### Result

Cytoplasm and Nucleolus	Green
Most areas of nucleus	Colorless

## Note

1. The blood smear or bone marrow smear should be uniform in thickness to avoid affecting the staining effect.
2. Blood cell smear staining requires fresh whole blood or EDTA anticoagulant.
3. After incubation with Acid Differentiation Solution, wash thoroughly, otherwise it will interfere with the staining of fast green.
4. If the dye color is too deep, it can be decolorized properly with methanol or alcohol, and it is better not to be re-dyed.
5. The pH value has certain influence on the dyeing. The slide should be clean and free of acid and alkali pollution to avoid affecting the dyeing effect.
6. For your safety and health, please wear experimental clothes and disposable gloves.