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Sperm Morphology Stain (Diff-Quick Method)

Cat: G2572

Size: 3×20mL/3×100mL

Storage: RT, avoid light, valid for 1 year.

Kit Components

Reagent	3×20mL	3×100mL	Storage
Reagent(A): Diff-Quick Fixative	20mL	100mL	RT
Reagent(B): Diff-Quick I	20mL	100mL	RT, avoid light
Reagent(C): Diff-Quick II	20mL	100mL	RT, avoid light

Introduction

Diff-Quick Staining is a fast staining method improved base on Wright Staining. It is one of the common staining methods in cytology test. The dye solution is made by the fast staining method recommended by WHO. Similar to Wright Stain, it is improved base on Romanowsky Stain technology principle. The staining result is very similar to that of Wright Stain, but Diff-Quick Stain takes a very short time, usually within 90 seconds. Diff Quick Stain can be used to evaluate sperm morphology.

Sperm Morphology Stain (Diff-Quick Method) is a sperm morphology staining solution recommended by WHO. Its dyeing principle is that sperm and proteins with different isoelectric points in cells have different charges under the same acidity. They can selectively combine with corresponding dyes to dye; the amino group dissociated by eosinophilic protein has positive charge, which can mix with acid dyes containing negative charge (such as eosinophilic xanthene and eosin) and be dyed red; the carboxyl group dissociated by basophil protein has negative charge, which can combine with basic dyes containing positive charge (basophil thioazabenzene and methylene blue) and be dyed blue; the number of amino group with positive charge and carboxyl group with negative charge dissociated by the neutral protein are the same. They simultaneously combine with acid dye and basic dye to be dyed purplish red, but the color is weak due to the equal dissociation charge. The dye is mainly used for sperm (cell) morphology evaluation, including fixative, which is very suitable for batch dyeing, and the background is clear without sediment. The dye is only used in scientific research, not for clinical diagnosis or other purposes.

Self Provided Materials

Slide, Distilled water, Microscope.

Protocol(for reference only)

1. Prepare smear: wash 2 slides thoroughly, then wash with 70% ethanol and dry. Add 5-20 μ L onto the slide, drag a drop of semen on the surface of the clean slide with the edge of the second slide, and make a smear. If the sperm density is too high, can dilute with physiological saline.
2. Put the smear in Diff-Quick Fixative or dry naturally, and fix for 15-20s.
3. Place the smear vertically on the absorbent paper to remove the excess liquid.
4. Stain with Diff-Quick I for 10-20s and place the smear vertically on the absorbent paper to remove the excess liquid.
5. Stain with Diff-Quick II for 5-10s and place the smear vertically on the absorbent paper to remove the excess liquid.
6. Wash with running water for 10-15 times to remove excess dye solution.
7. Place the smear vertically on the absorbent paper to remove the excess liquid and make it completely dry. View under the microscope.

Result

Acrosome region of sperm head	Light Blue
Postacrosomal region	Deep Blue
Middle part	Maybe Slight Red
Rear part	Blue or Slight Red

Note

1. If the sperm density is more than 20×10^6 / mL, take 5 μ L; if the sperm density is less than 20×10^6 / mL, take

5-20 μ L.

2. Diff-Quick Fixative is a special fixative for Sperm Diff-Quick Staining Method. If use a lot ,can replace with methanol.
3. The sperm sample should be fresh, and the sperm smear should be even and thick, so as not to affect the staining effect. If the sperm density is too high, adjust the sperm density with physiological saline.
4. Avoid removing the dye solution directly or wash the smear vigorously. Do not discard the dye first, so as to prevent the dye from settling on the smear.
5. The dye solution can be reused, but it can not be repeated many times. If there is sediment, filter before use.
6. If the dyed color is too deep, it can be decolorized properly with methanol or alcohol, and it is better not to be re-dyed.
7. If the dyed color is too deep or too shallow, adjust the dyeing time or working solution concentration.
8. The pH value has certain influence on the dyeing. The slide should be clean and free of acid and alkali pollution to avoid affecting the dyeing effect.
9. For your safety and health, please wear experimental clothes and disposable gloves.