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规格: 2×100mL/2×500mL

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Cat: G3100

Size: 2×100mL/2×500mL

Storage: RT, avoid light, valid for 2 years.

Kit Components

Reagent	2×100mL	2×500mL	Storage
Reagent(A): May-Grunwald Stain	100mL	500mL	RT, avoid light
Reagent(B): MG Phosphate Buffer	100mL	500mL	RT
Mix Reagent A and B in an equal amount to form May-Grunwald Working Solution. It is not suitable to prepare in advance.			

Introduction

May-Grunwald Stain is often used for staining of tissue sections, blood and cell smears, bacteria, chromosome banding, protozoa parasites and so on. It is especially for helicobacter pyloric. Helicobacter pylori is also known as campylobacter pylori. It has been proved that the bacteria is closely related to chronic gastritis and peptic ulcer. Helicobacter pyloric in the stomach is generally arc-shaped, S-shaped or seagull shaped, and sometimes 3-4 bends are seen in spiral shape, often in the shape of fish. Helicobacter pyloric is mostly found between the surface epithelium and the mucosal layer of the gastric mucosa, and close to the surface epithelial cells, some of them enter into the cytoplasm of the epithelial cells, and the bacteria are also found in the gastric fovea and the superficial glandular cavity of the mucosa.

The main staining methods of helicobacter pyloric are Methylene Blue method, Silver Nitrate method, May-Grunwald-Giemsa (MGG) method and Basic Fuchsin method. Silver Nitrate method has a clear contrast and the stained section can be preserved for a long time, but the operation is more time-consuming and troublesome. Other methods are simple and convenient, but the stained section is easy to fade. After May-Grunwald staining, helicobacter pyloric is blue, collagen fiber is red, red blood cell is green, gastric mucosa epithelium is light blue, the nucleus is dark blue.

Self Provided Materials

10% formalin fixative, distilled water, absolute ethanol

Protocols(for reference only)

1. Fix the tissue in 10% formalin and conventionally dehydrate and embed.
2. Cut into sections in 4μm thick, and dewax to distilled water.
3. Wash with distilled water once.
4. Wipe the water around the section, drop May-Grunwald Working Solution onto the section and stain for 10-15min, and then discard the solution.
5. Quickly wash away the remaining solution with absolute alcohol and dry it slightly.
6. Transparent by xylene and seal with resinene.

Result

Helicobacter pyloric, Nucleus, Epithelium of gastric mucosa	Blue
Collagen fiber	Red
Red blood cell	Green

Note

1. May-Grunwald Working Solution must be prepared before use and it is not easy to preserve.
2. Wash quickly in the final absolute ethanol, otherwise it is easy to decolorize.
3. Helicobacter pyloric can be showed when only staining with Giemsa Stain and the color is clear. But the best effect is staining with combination of May-Grunwald Stain and Giemsa Stain.
4. For your safety and health, please wear experimental clothes and disposable gloves.