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## Phospholipid-Fe Hematoxylin Stain Solution

**Cat:** G3250

**Size:** 3×100mL

**Storage:** RT, avoid light, valid for 6 months.

### Kit Components

Reagent		3×100mL	Storage
Reagent(A): FeH Fixative		100mL	RT, avoid light
Reagent(B): FeH Stain Working Solution	Reagent(B1):FeH Hematoxylin	5×0.1g	RT, avoid light
	Reagent(B2):Hematoxylin Solvent	25mL	RT
	Reagent(B3):FeH Buffer	75mL	RT, avoid light
Reagent(C): FeH Differentiation		100mL	RT

### Introduction

Phospholipid is a kind of compound lipid containing phosphoric acid. Phospholipids are the components of biofilms, which are divided into glycerophospholipids and sphingolipids. Phospholipids are amphoteric molecules. One end is hydrophilic nitrogen or phosphorus containing head, and the other end is hydrophobic long alkyl chain. Elleder found that phospholipids can be displayed by ferritin, which is more simple, rapid and sensitive than Dah. However, acetone should be used to degrease the sample properly, and the phospholipid and nucleus will be stained at the same time.

### Self Provided Materials

Chloroform, Acetone, Methanol, Distilled Water

### Protocol(for reference only)

1. One slice was soaked in chloroform methanol solution (1:1) for 1 h as negative control, and another slice was soaked in acetone precooled at 2-8 °C for 15 min.
2. The rats were fixed with FEH for 30 min. Rinse with distilled water twice, 3 min each time.
3. To prepare FEH Stein in advance, add 0.1g of FeH Hematoxylin to 5ml of Hematoxylin Solvent. After fully dissolving, add 15ml FeH Buffer to form the FeH Stain Working Solution.
4. Staining with FEH stain for 7-10 min. Rinse with distilled water twice, 30s each time.
5. After the FEH differentiation solution was diluted with distilled water, the slices were immersed in the diluted FEH differentiation solution for several times. Rinse with tap water.
6. Dehydrated with 95% ethanol, transparent with xylene and sealed with neutral resin.

### Result

Phospholipid	Blue
Nucleus	Blue

### Note

1. Phospholipids are easy to dissolve, so the tissue should be fixed immediately after removal, otherwise it is difficult to stain.
2. The results showed that the staining solution had good staining effect on frozen section.
3. The stained specimen must be kept away from light, otherwise it is easy to fade.
4. For your safety and health, please wear lab clothes and disposable gloves.