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Sodium Sulfate Solution, 5%, Water Solvent

Cat: G3420

Size: 500mL

Storage: RT, valid at least for 6 months.

Introduction

In the process of tissue sectioning, when some tissues contain bone or calcification, it is not suitable to directly use paraffin embedded sections for calcium containing tissues. This is because the density of calcium and paraffin is different, so it is difficult to cut a complete section. It is better to fix the calcified tissue before decalcification or both. In the process of decalcification, acids break the cross chain of protein molecules and release hydrophilic groups which can combine with water molecules, causing tissue expansion. After decalcification, the acid in the tissue should be neutralized immediately to make the tissue deacidification and neutralization quickly, so as to ensure the quality of the section and the staining effect.

Sodium Sulfate Solution, 5%, Water Solvent is a commonly used post-treatment solution for bone tissue decalcification, which is usually recommended to be treated overnight.

Protocol(*for reference only*)

1. When decalcifying the bone tissue, the thickness of the material should not be too thick, usually no more than 5mm.
2. After decalcification with acidic decalcification solution, the tissue was soaked in distilled water for three times, 1 min each time.
3. The tissue was transferred to Sodium Sulfate Solution, 5%, Water Solvent for 24 h.
4. Wash twice with distilled water, 1h each time, transfer to 70% ethanol for dehydration and embedding.

Note

1. Proper heating can speed up the treatment, but high temperature will also damage the microstructure, so it is recommended to control the temperature below 40 °C.
2. Long time washing may cause tissue swelling, so it is recommended to transfer it to 70% ethanol for preservation or dehydration embedding within 6h.