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规格: 3×30mL/3×100mL/3×250mL/3×500mL

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- [1]YanxiaGao.HDAC1 promotes artery injury through activation of VAV3 by binding to miR-182-5p in atherosclerotic mice model.CeLLuLar SignaLLing.November 2020.(IF 3.968)
- [2]Xifeng Yang Yao Sun Ying Zhang Shan Han,et al.Downregulation of miR- 181b inhibits human colon cancer cell proliferation by targeting CYLD and inhibiting the NF- κB signaling pathway.International Journal of Molecular Medicine.September 2020.(IF 3.098)

## Quick Giemsa Stain Kit

**Cat:** G4640

**Size:** 3×30mL/3×100mL/3×250mL/3×500mL

**Storage:** RT, avoid light, valid for 1 year.

### Kit Components

Reagent	3×30mL	3×100mL	3×250mL	3×500mL	Storage
Reagent(A): Giemsa Stain Solution	30mL	100mL	250mL	500mL	RT
Reagent(B): Giemsa Buffer	2×30mL	2×100mL	2×250mL	2×500mL	RT

### Introduction

Giemsa dye is composed of azure and eosin. all kinds of components present different colors and can be distinguished.

The rapid Giemsa staining solution is made of imported Giemsa pigment dye. Without pre fixation, it can dye the nucleus into purplish red or blue purple, and the cytoplasm into pink, showing clear cell staining images under light microscope. It is mainly used to show the difference of various blood cell morphology and size in blood smear, with good staining effect, strong staining power and clear staining.

### Protocol (for reference only)

1. Prepare the blood smear in usual method. After the smear dried, place on the dyeing rack, and add 200-400μl Giemsa stain solution to cover the blood smear quickly. Staying for 90-120s to fix and pre-dye the smear.
2. Add 1-2 times the volume of Giemsa buffer directly, shake the slices lightly or blow air with ear ball to make the dye well mixed, and stain for 5-8 min.
3. Slowly rinse the glass slide with tap water for 30s (pay attention not to pour out the dye solution or directly rinse the blood film), dry it and then examine it under microscope.

### Result

The Red Blood Cells	light red, and the center is slightly pale and slightly discoid.
The Lymphocyte	Membrane is blue purple, the nucleus is purplish red, the cytoplasm was light red, and the granules were clear.

### Note

1. For your safety and health, please wear lab clothes and disposable gloves.
2. Smear thickness should be moderate, too thick or too thin is not easy to get the result wanted.
3. Giemsa Stain Solution is volatile and should not be too little. It is recommended to add 300-400ul for each slice. The reagent is prepared by conventional method and filtered by filter paper, so as to avoid the failure of the reaction between the dye solution and water.
4. PH had effect on cell staining. The slides used for staining must be clean and free from acid and alkali pollution, so as not to affect the staining results.

### Reference

- [1] Yanxia Gao. HDAC1 promotes artery injury through activation of VAV3 by binding to miR-182-5p in atherosclerotic mice model. *Cellular Signaling*. November 2020. (IF 3.968)
- [2] Xifeng Yang Yao Sun Ying Zhang Shan Han, et al. Downregulation of miR-181b inhibits human colon cancer cell proliferation by targeting CYLD and inhibiting the NF-κB signaling pathway. *International Journal of Molecular Medicine*. September 2020. (IF 3.098)