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**Cat:** G1135

**Size:** 3×50mL

**Storage:** 2-8℃, avoid light, valid for 6 months.

### Kit Components

Reagent	3×50mL	Storage
Reagent(A): Leifson buffer	50mL	2-8℃
Reagent(B): Leifson Diluent	50mL	2-8℃
Reagent(C): Fuchsin Solution	50mL	2-8℃, avoid light
Before use, mix Reagent (A):(B):(C) with the ratio of 1:1:1 to form Leifson Staining Solution. Use it immediately, not for a long time.		

### Introduction

Bacterial flagellum is the motor organ of bacteria. *Helicobacter pylori* can pass through the mucus layer of gastric epithelial cells from the strongly acidic gastric cavity to the neutral environment of gastric epithelial cells, which is a good example of flagellum movement. Flagella morphology, quantity and location of flagella in bacteria can be observed by flagella staining. Flagellum quantity and its distribution on bacteria are one of the important bases for identifying bacteria.

This kit uses Leifson method, also known as Leifson's flagella stain method. The reagent is more sensitive, the operation is simple and the result is more reliable.

### Self Provided Materials

Alcohol Lamp, Slide, Distilled Water, Inoculation Ring, Microscope

### Protocol(for reference only)

- Drop 2 drops of distilled water on a clean, fat-free slide.
- Select the sterile distilled water with the inoculation ring, then contact the colony on the blood plate, allow the bacteria to swim into the inoculation ring distilled water, and then move the inoculation ring to the top of the distilled water on the slide twice.
- Shake the slide to distribute the bacteria evenly. Do not grind and stir with the slide to avoid flagella falling off.
- Drip Leifson Staining Solution on the slide. When more than half of the area of the Leifson Staining Solution appears metallic gloss film, generally dye for 3-10 min. Gently wash the staining solution with distilled water and air-dry naturally.
- Microscopic examination should begin at the edge of the smear and gradually move to the center from the outside to the inside. Where bacteria are less distributed, flagellum is easy to observe. Where bacteria are concentrated, flagellum is blocked by fungi and is not easy to observe.

### Result

Thalli and Flagellum	Red in Different Degree
----------------------	-------------------------

### Note

- Slides should be clean and free from oil pollution.
- The bacteria sample for staining must be in the active stage of growth, which should be successive subculture for many times.
- Operate carefully in the dyeing process to prevent flagella falling off.
- The prepared Leifson Staining Solution should be used within 2 h as far as possible avoid storing for a long time.
- It is not appropriate to fix with high-temperature flame.
- For your safety and health, please wear experimental clothes and disposable gloves.