

卢戈碘液（标准碘液）

货号：G1069

规格：100mL/500mL

保存：2-8°C，避光保存，有效期 3 个月。

产品介绍：

卢戈碘液（标准碘液）又叫鲁哥氏染色液(Lugol's 碘液,5%)，是碘和碘化钾的水溶液。5 克碘(I₂)和 10 克碘化钾(KI)溶于 85mL 蒸馏水中。常应用于检测叶片、食物等样品中的淀粉，可将淀粉染成黑色或蓝黑色。用于细胞染色，使细胞核着色从而更易观察。

相关产品：

G1066 卢戈碘液（革兰氏染色）

G1062 结晶紫染色液（1%）

G1063 结晶紫染色液（0.1%）

G1060 革兰氏染色液(试剂盒)

Lugol's Iodine Solution(Iodine Standard Solution)

Cat: G1069

Size: 100mL/500mL

Storage: 2-8°C, avoid light, valid for 3 months.

Introduction

Lugol's Iodine Solution(Iodine Standard Solution) also known as Lugol's dye (5%), is an aqueous solution of iodine and potassium iodide. 5 grams of iodine (I₂) and 10 grams of potassium iodide (KI) were dissolved in 85mL distilled water. It is often used in the detection of starch in leaves, food and other samples. The starch can be dyed black or blue-black. If used for cell staining, please stain the nucleus to make it easier to view.

Related Products

G1066 Lugol Solution(Gram Method)

G1062 Crystal Violet Stain Solution, 1%

G1063 Crystal Violet Stain solution, 0.1%

G1060 Gram Stain Kit