

鞣酸水溶液(5%)

货号: G3440

规格: 100mL

保存: 室温, 避光保存, 有效期 1 年。

产品介绍:

鞣酸又称单宁酸(Tannic acid), 多用于鞣革制剂, CAS 号为 1401-55-4, 分子量为 1701.20。

鞣酸水溶液(5%)由鞣酸、去离子水组成, 具有媒染效果, 常与碱性染料联用。可用于肌纤维染色、鞭毛染色和半薄切片染色和微血管染色。

操作步骤: (仅供参考)

1. 根据实验具体要求操作。

注意事项:

1. 试剂均有刺激性, 请小心操作。
2. 为了您的安全和健康, 请穿实验服并戴一次性手套操作。

Tannic Acid Solution, 5%, Water Solvent

Cat: G3440

Size: 100mL

Storage: RT, avoid light, valid for 1 year.

Introduction

Tannic acid, also known as tannic acid, is mainly used in tanning preparations. The CAS number is 1401-55-4, and the molecular weight is 1701.20.

Tannic Acid Solution, 5%, Water Solvent is composed of tannic acid and deionized water. It has mordant effect and is often combined with basic dyes. It can be used for muscle fiber staining, flagella staining, semi thin section staining and microvascular staining.

Protocol (for reference only)

Operate according to the specific requirements of the experiment.

Note

1. The solution is irritant, please operate carefully.
2. For your safety and health, please wear experimental clothes and disposable gloves.