

## Mouse Anti Rabbit IgG Monoclonal Antibody

Cat. Number: K20003M

Quantity size: 1mg/5mg/10mg

**Concentration:** 2mg/ml buffer=0.01M PBS (pH7.2-7.4)with 0.03% Proclin300.

**Storage:** Shipped at 4°C, Store at -20°C (Avoid repeated freeze/thaw cycles).

## **Background:**

Immunoglobulin G (IgG), is one of the most abundant proteins in serum with normal levels between 8-17 mg/mL in adult blood. IgG is important for our defence against microorganisms and the molecules are produced by B lymphocytes as a part of our adaptive immune response. The IgG molecule has two separate functions; to bind to the pathogen that elicited the response and to recruit other cells and molecules to destroy the antigen. The variability of the IgG pool is generated by somatic recombination and the number of specificities in an individual at a given time point is estimated to be 1011 variants.

**Specificity:** Mouse monoclonal IgG.

Species reactivity: Rabbit.

Immunogen: Rabbit IgG full length protein.

Purification: Affinity purification.

Predicted Molecular Weight: 150kDa.

Clone ID:6C6

## **Application:**

WB 1:1:1000-20000. ELISA:1:5000-50000.

Not yet tested in other applications.

Optimal working dilutions must be determined by the end user.

## **Important Note:**

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.