

Anti-SPIB Polyclonal Antibody

Cat: K109500P


Summary:

【Product name】 : Anti-SPIB antibody	【Source】 : Rabbit
【Isotype】 : IgG	【Species reactivity】 : Human Mouse Rat
【Swiss Prot】 : Q01892	【Gene ID】 : 6689
【Calculated】 : MW:19/30kDa	【Observed】 : MW:29kDa
【Purification】 : Affinity purification	
【Tested applications】 : WB IHC	
【Recommended dilution】 : WB 1:2000-5000. IHC 1:25-100.	
【WB Positive sample】 : PC-3M-IE8, THP-1	
【IHC Positive sample】 : Human stomach cancer	
【Subcellular location】 : Cytoplasm	
【Immunogen】 : Recombinant protein of human SPIB	
【Storage】 : Shipped at 4°C. Upon delivery aliquot and store at -20°C	


Background:

The protein encoded by this gene is a transcriptional activator that binds to the PU-box (5'-GAGGAA-3') and acts as a lymphoid-specific enhancer. Four transcript variants encoding different isoforms have been found for this gene.

Verified picture



Western blot analysis with SPIB antibody diluted at 1:4000; Lane: PC-3M-IE8, THP-1



Immunohistochemistry of paraffin-embedded Human stomach cancer with SPIB antibody diluted at 1:40